

# Investigating Oxidation- reduction Reactions

## Chemistry 101: *General Chemistry*

Post-Lab & Lab Report #7



*Sitthiphol Yuwanaboon*

*Professor Nina Ram*

*Lab Professor*

*January 29, 2020*

## Investigating Oxidation- Reduction Reactions

**Purposes:** To investigate the chemical reaction (Redox) for reactants which is oxidize and Reduced for looking chemical forming such precipitation, gases ,and color changes. This is type of things we are looking for to define what happen that cause the reactants become oxidize and reduced. That means the electron has lost and gained in the chemical reaction. For part IV and V to find the amount of chemical element in the in the reaction which is reacted.

### Procedure:

**Part I:** calculate the mass of  $\text{KMnO}_4$  to make an .02 M solution of  $\text{KMnO}_4$  dilute with the water. The solution require a mass of .94824 gram of potassium permanganate then dissolve in boiling water about 300 mL of the solution we should get for the experiment in the next part. Then save it in the container of 1000 mL

### Part III:

Calculate the mass require of oxalic acid to react with  $\text{KMnO}_4$  in process of titration to find the concentration. Put in the flask and add the 50 mL then and add the sulfuric acid to the solution about 13-14 mL. heat with hot plate till it reach 80-90 Celsius then titrate with  $\text{KMnO}_4$  the solution you make for the experiment. Wait for the color change to light pinkish then you stop titrate then you measure the quantity of  $\text{KMnO}_4$  you added then calculate the morality of solution. Repeat for two other flasks.

### Part IV:

Instead of oxalic acid we put some unknown sample for the instructor. We can use the same process as above weight out sample and divide in four parts then repeat Step from Part III

### Part V:

Instead of oxalic acid we put some unknown sample for the instructor. We can use the same process as above weight out sample and divide in four parts then repeat Step from Part III to find out the ferrous ion in the solution and percentage of by mass used in the experiment and how much are there. Steps are same in part III but it is just a new unknown which we got from the professor.

**Data charts:**

Procedures	Observations	Conclusions
<b>Test for sulfate</b> 1. Dissolve small $\text{Na}_2\text{SO}_4$ in 2 mL deionized water in the test tube. Then put 20 drop of 6M $\text{HCl}$ . No more than 4-6 drop of $\text{BaCl}_2$	There is precipitate formed in test tube like gelatinous white matter on the top of solution	$\text{SO}_4^{2-}$ is present because there is precipitate in the solution.
<b>Behavior of <math>\text{HSO}_3^-</math></b> 1. 2 ml $\text{NaHSO}_3$ then test for sulfate but add 2 ml distilled water to bisulfite solution watch out the precipitate in test tube	No precipitate in the solution color is transparent. No color	There is no chemical reaction occur because there redox of $\text{SO}_4^{2-}$ and no precipitate in the test tube
<b>A Effect of <math>\text{HSO}_3^-</math> on <math>\text{MnO}_4^-</math></b> 1. 2 ml $\text{NaHSO}_3$ then put the dilute $\text{HCl}$ then add a $\text{KMnO}_4$ one at a time and mixing every drop and see color changes then test for sulfate	The purple color changed to clear solution  Cloudy are form when we add $\text{BaCl}_2$	There is $\text{Mn}^{2+}$ pink and it is too dilute to see and redox and $\text{SO}_4^{2-}$ is present because of precipitate in solution
<b>B:</b> 1. 2 ml $\text{NaHSO}_3$ then put 1% 2-3 of Hydrogen peroxide then test for sulfate.	Cloudy less ,little precipitate	Little reaction , less reactant , little precipitate or none. Little chemical reaction
<b>C:</b> 1. 2 ml $\text{NaHSO}_3$ then put $\text{Br}_2$ 14-15 drop see color changes then test for Sulfate	Tiny cloudy and little precipitate the color $\text{Br}_2$ is yellow	Little reaction , less reactant , little precipitate or none. Little chemical reaction
<b>D:</b> 1. 2 ml $\text{NaHSO}_3$ then put dilute $\text{HCl}$ then add $\text{K}_2\text{Cr}_2\text{O}_7$ solution one drop at a time mixed after each drop then test for sulfate	No color changed when we add $\text{K}_2\text{Cr}_2\text{O}_7$ in solution Nothing seem to happen. No color change when test with sulfate	No chemical reaction occur because there is no precipitate are formed when we test with $\text{SO}_4^{2-}$

Procedures	Observations	Conclusions
<b>3a:</b> 1. to first tube add 5 drop dilute $\text{H}_2\text{SO}_4$ then put Hydrogen peroxide 8-10 drop see the bubble are formed	Don't see any gas, by condition we have	There is decomposition chemical reaction because it form a gas
<b>3B:</b> 1. to first tube add 5 drop dilute $\text{H}_2\text{SO}_4$ then add ammonium sulfate. See color change.	Color become clear from the purple color.	There is ion transform because its color change.
<b>3C:</b> 1. to first tube add 5 drop dilute $\text{H}_2\text{SO}_4$ . add small amount of KI then see color change. Then do add sodium thiosulfate until color disappear. See what color has changed	When add KI the color become yellow dark similar to test tube I2  The color change when we add the starch on it to brownish then we add $\text{Na}_2\text{S}_2\text{O}_4$ the test tube become clear again	I2 which make the test tube look yellow. The chemical reaction occur it color change. Then acid change it color to restate of phase.
<b>2B:</b> 1. 2 ml $\text{NaHSO}_3$ then put 5-6 drop dilute $\text{NaOH}$ and 2 drop of $\text{KMnO}_4$ solution see color change	Light yellow formed when we put $\text{KMnO}_4$ to the solution Look like more yellow.	There ion which change the color to yellow
<b>2C:</b> 1. 2 ml $\text{NaHSO}_3$ then put $\text{Na}_2\text{CO}_3$ stir gently till it dissolved then add 2 drop of $\text{KMnO}_4$	For a second we can see yellow brown color we add $\text{KMnO}_4$	Little chemical reaction occur because we don't see big amount of changes .

Procedures	Observations	Conclusions
<b>5:</b> 1. Dissolve oxalic acid and then add 20 drop dilute of $\text{H}_2\text{SO}_4$ . Separate in two test tube one with hot burner and at room temperature. After 5 minute then add 3-4 drop $\text{KMnO}_4$ to each test tube which goes faster in reaction.	With heat dissolve faster than no heat applied. The become transparent we put the $\text{KMnO}_4$ but the color of cold one is still purple on the of it.	Temperature cause chemical reaction to change phase faster and go back in the same stage( Shift)

### Part III: Standardization of $\text{KMnO}_4$

	<b>Trial 1</b>	<b>Trial 1</b>	<b>Trial 3</b>
<b>Mass of vial</b>	6.1928g	6.2471g	6.1834g
<b>Mass of vial with acid</b>	6.2871g	6.3420g	6.2778g
<b>Mass after transfer</b>	6.1928g	6.2475g	6.1840g
<b>Mass of oxalic acid</b>	.0943g	.0945g	.0938g
<b><math>\text{KMnO}_4</math> Initial buret reading</b>	0.0	12.0	24.1
<b><math>\text{KMnO}_4</math> final buret reading</b>	12.0	24.0	36.1
<b>Total Volume added</b>	12.0	12.0	12.0
<b>Mole of Oxalic acid</b>	$7.479 \times 10^{-4} \text{ mol}$	$7.4952 \times 10^{-4} \text{ mol}$	$7.4397 \times 10^{-4} \text{ mol}$
<b>Mole of <math>\text{KMnO}_4</math></b>	$2.992 \times 10^{-4} \text{ mol}$	$2.9998 \times 10^{-4} \text{ mol}$	$2.975 \times 10^{-4} \text{ mol}$
<b><math>\text{KMnO}_4</math> concentration</b>	.024933 mol/ L	0.024998 mol/ L	0.0247916 mol/ L
<b>Average concentration</b>	.024907 mol/ L	.024907 mol/ L	.024907 mol/ L
<b>Sig Average concentration</b>	0.0249 mol / L	0.0249 mol / L	0.0249 mol / L
<b>Deviation from average</b>	0.000026	0.000091	0.000115
<b>Relative Deviation</b>	1.04	3.65	4.63
<b>Average relative deviations</b>	3.11	3.11	3.11

From the beginning, we measure masses accurately and precision of titration is also accurate but the dilute solution significant is high from the reference number that we can refer to. Relative deviation and the different relative deviations from average are less than 10 ppt. so we measure and repeat the experiment accurately and there is might be something wrong with concentration of the  $\text{KMnO}_4$  dilute solution. If I get around 15 mL of  $\text{KMnO}_4$  the concentration would be accurate for what we expected from the experiment.

#### Part IV: Percent Oxalate Ion

	<b>Trial 1</b>	<b>Trial 1</b>	<b>Trial 3</b>
<b>Mass of vial</b>	6.1860g	6.2510g	6.1849g
<b>Mass of vial with unknown</b>	6.4480g	6.5128g	6.4455g
<b>Mass after transfer</b>	6.1860g	6.2524g	6.1849g
<b>Mass unknown used</b>	.2620g	.2618g	.2606g
<b>KMnO<sub>4</sub> Initial buret reading</b>	0.0	0.0	0.0
<b>KMnO<sub>4</sub> final buret reading</b>	28.7	28.5	28.5
<b>Total Volume added</b>	28.7	28.5	28.5
<b>Mole of Oxalate Ion</b>	$1.7871 \times 10^{-3} \text{ mol}$	$1.7745 \times 10^{-3} \text{ mol}$	$1.7745 \times 10^{-3} \text{ mol}$
<b>Mole of KMnO<sub>4</sub></b>	$7.1483 \times 10^{-4} \text{ mol}$	$7.098 \times 10^{-4} \text{ mol}$	$7.098 \times 10^{-4} \text{ mol}$
<b>Mass of Oxalate Ion</b>	.15729g	.15619g	.15619g
<b>% Percent Oxalate Ion</b>	<u>60.03</u>	<u>59.66</u>	<u>59.94</u>
<b>Average % Percent Oxalate Ion</b>	59.87763%	59.87763%	59.87763%
<b>Sig Average % Percent Oxalate Ion</b>	59.9%	59.9%	59.9%
<b>Deviation from Average</b>	.1523	.2176	.05236
<b>Relative Deviation</b>	2.543	3.634	1.040
<b>Average Relative Deviations</b>	2.405	2.405	2.405
<b>Sig Average Relative Deviations</b>	2.41	2.41	2.41

I would not round much of for significant because the Relative deviation will get larger if I round up the number. So I put the sig Avg for the percent oxalate in the solution of the experiment. So we make it better that way. And accuracy and precision are in the right spot. PPT is less than 10, so precision is precise and accurate. So there will be only three sig fig in experiment, that the end of part IV by find the mass of oxalate but we have to find the mole ratio of **KMnO<sub>4</sub> and Oxalate is 2:5** so we time 5 divide by two, so we will the number of mole of oxalate then convert gram of oxalate. Then find the percent ratio in trials we got and weigh it.

### Part V: Percent Ferrous Ion in Redox Experiment

	<b>Trial 1</b>	<b>Trial 1</b>	<b>Trial 3</b>
<b>Mass of vial</b>	6.1859g	6.1859g	6.1859g
<b>Mass unknown used</b>	1.0048g	1.0053g	1.0059g
<b>KMnO<sub>4</sub> Initial buret reading</b>	0.0	15.9	22.6
<b>KMnO<sub>4</sub> final buret reading</b>	15.9	31.8	38.5
<b>Total Volume added</b>	15.9	15.9	15.9
<b>Mole of unknown Ion</b>	$1.59 \times 10^{-3}$ mol	$1.59 \times 10^{-3}$ mol	$1.59 \times 10^{-3}$ mol
<b>Mole of KMnO<sub>4</sub></b>	$3.18 \times 10^{-4}$ mol	$3.18 \times 10^{-4}$ mol	$3.18 \times 10^{-4}$ mol
<b>Mass of Unknown</b>	0.0888015g	0.0888015g	0.0888015g
<b>% Percent ferrous Ion</b>	8.838%	8.833%	8.828%
<b>Average % Percent Ferrous Ion</b>	8.833%	8.833%	8.833%
<b>Sig Average % Percent Ferrous Ion</b>	8.83%	8.83%	8.83%
<b>Deviation from Average</b>	.005	.0	.005
<b>Relative Deviation</b>	.5661	0	.5661
<b>Average Relative Deviations</b>	.3773	.3773	.3773
<b>Sig Average Relative Deviations</b>	$3.77 \times 10^{-1}$	$3.77 \times 10^{-1}$	$3.77 \times 10^{-1}$

We divide a mass of 4.0130 with a vial four parts then we got the volume added then covert to mole using the morality of Solution we got from part III then time by mole ratio of ferrous ion 5:2 of **KMnO<sub>4</sub>** so we time 5 divide 2 then we get moles of unknown then using molar mass of Fe :55.85 g/ mol to convert to gram of Ferrous ion then we compare the percentage of what we have from calculation to the what we used in the solution of mass unknown per trial.



**Conclusion:**

	<b>Average value</b>	<b>Average relative deviation</b>
<b>Part III: Average molarity of KMnO<sub>4</sub></b>	.0249 mol/ L	3.11
<b>Part IV: unknown number</b>	1039R (unsure)	
<b>Sig Average % Percent Oxalate Ion</b>	59.9%	2.41
<b>Part V:unknown number</b>	760R (unsure)	
<b>Sig Average % Percent Ferrous Ion</b>	8.83%	$3.77 \times 10^{-1}$